

Flushing Kit

For All Pumps

The Masport Flushing Kit has been specifically designed by Masport for the maintenance of all Masport Pumps.

The Masport Flushing Kit provides a simple solution to be able to perform regular flushing of your Masport Pump. Flushing the pump on a regular basis helps to extend the life of the pump and ultimately save you money in the long run.

Masport quality and performance

- ▶ **Extends pump life** - keep pumps clean and running for years
- ▶ **Durable design** - UV resistance reservoir and heavy duty fittings
- ▶ **Easy to install** - fittings provided for quick installation on all Masport pumps
- ▶ **User friendly** - detailed operating instructions on every flushing reservoir
- ▶ **Steel and PVC reservoirs available**
- ▶ **Flushing Fluid specially formulated** for Masport Pumps available
- ▶ **1 Gallon Flushing Fluid Bottles**
Part # 13989
- ▶ **6 Gallon Flushing Fluid Cases**
Part # 13990



Flushing Kit


For All Pumps



Flushing Kit

Instructions

Exploded View



No.	Description	Code	Qty
1	Tank Cap	36002	1
2	Lubricator Tank	36001	1
3	Hose 1/4"	24299	2
4	Hose Clamp	24310	1
5	Pushlock 1/4"	27001	3
6	Elbow 1/4"	27002	2
7	Flow Restrictor	27032	1
8	Bracket	33461	1
9	Nipple 1/4"	30000	1
10	Ball Valve	25154	1

Flushing Procedure

Flush weekly, preferably at the end of the work week, or anytime carry-over contaminates the pump.

- ▶ Start vacuum pump and bring vacuum level to 15" Hg.
- ▶ Open ball valve to begin flushing fluid flow.
- ▶ Continue to run pump until fluid level drops one graduation (approximately 7 fl oz).
- ▶ Close ball valve
- ▶ Continue to operate pump one additional minute to ensure all fluid is purged from pump.
- ▶ Stop vacuum pump.
- ▶ Drain Oil Separator and dispose used fluid at approved recycling site.

Never reuse flushing fluid or pump damage may occur.

